



This is a numbered callout. It is not part of the page design, instead it typically refers to a comment, question, or function.

ReadyAdminConfig.100

This is the file name of the mockup (for reference purposes only).




This is a yellow sticky note. It is not part of the page design. It includes comments. Note: If pink, it is a question for the customer.



This is a pink sticky note. It is not part of the page design. It is a question for the customer.

Note: the web page URL is also used to indicate what page you are on.

Site Map






READY Measurement Administration

Username:

Password:

A Web Page

http://ReadyAdmin/site-listing

READY Measurement Administration [Logout](#)

Site Management Configuration Management

Name	Password	Action
osu	osu_password	Delete
uab	uab_password	Delete

[Create Site Account](#)

ReadyAdminSite.100

- 1 Listing of site accounts
- 2 To edit, select name (hyperlink)
- 3 To create a new account for a site, select "Create Site Account"

For the most part, this is the new detail page for a site (i.e., some fields are filled in to demonstrate functionality).

Is the Greeting site dependent or configuration dependent? If possible, would like one for both site and configuration (11/13/2014). This is a placeholder may or may not be in December deliverable as it is low priority

- 1 Dropdown menu that contains all available configurations in application.
- 2 Select a configuration and press the Add button to add that configuration to the configuration listing for the given account.
- 3 Select the default configuration for the given account.
- 4 Press the Delete button to remove that configuration from the listing for the given account.

- 5 Predictive search with autocomplete functionality. Providers obtained via a web service to the national provider's index. Once a provider is selected, they are added to the list of providers for the given account.
- 6 List of providers for the given account.
- 7 Press the Delete button to remove that provider from the listing for the given account.

Account

8

Configurations

2

Providers

Account

Save

Press the Delete button to remove that provider from the listing for the given account.

A Web Page

http://ReadyAdmin/configuration-listing

READY Measurement Administration [Logout](#)

Site Management Configuration Management

Name	Action
Configuration 1	Delete
Configuration 2	Delete
Configuration 3	Delete
Configuration 4	Delete
Configuration 5	Delete
Configuration 6	Delete
Configuration 7	Delete

Create Configuration

ReadyAdminConfig100

- 1 Listing of available configurations (all institutions).
- 2 To edit, select name (hyperlink)
- 3 To remove a configuration, press Delete.
- 4 To create a new Configuration, select "Create Configuration"

Configuration 2

☐ No


5

If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Sad faces to show in pictograph 3. Note: The number of neutral faces is $1000 - \#Sad$.

Thank you for accessing the Rheumatoid Arthritis Disease activity (READY) measurement tool. This electronic application can be used by both patients and physicians to assess how your RA disease is doing.

If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Sad faces to show in pictograph 3. Note: The number of neutral faces is 1000 - #Sad.

Configuration Diagnoses and Instruments



A screenshot of a web browser toolbar. The toolbar contains icons for Bold (B), Italic (I), Underline (U), Text Color (color swatch), Style (dropdown menu), List (bulleted list), Link (chain link), Unlink (broken chain link), Image (picture icon), and Smiley (smiley face icon). The 'Style' dropdown menu is open, showing a list of styles.

5 If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Sad faces to show in pictograph 3. Note: The number of neutral faces is $1000 - \#Sad$.

2

Name: Configuration 3

1

Save Results? ☐ Yes ☒ No

Configuration Diagnoses and Instruments

[illegible]

ReadyAdminConfig.107

This is the detail page for updating an existing configuration. It appears when the configuration is selected from the configuration listing).

This page takes advantage of AJAX calls, so updates are immediate (extra buttons for save and cancel are not necessary).

If Save Results is No, iPad app will NOT need Patient ID and only Not Sure checkboxes are enabled.

3

If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Happy faces to show in pictograph 1. Note: The number of neutral faces is $100 - \# \text{Happy}$.

4

If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Sad faces to show in pictograph 2. Note: The number of neutral faces is $100 - \#Sad$.

5

If checkbox in row above is selected, text box is enabled and numbers can be added for the number of Sad faces to show in pictograph 3. Note: The number of neutral faces is 1000 - #Sad.

Configuration Greeting

Thank you for accessing the Rheumatoid Arthritis Disease activity (READY) measurement tool. This electronic application can be used by both patients and physicians to assess how your RA disease is doing.

<http://ReadyAdmin/configuration-listing>

READY Measurement Administration

[Logout](#)

Site Management

Configuration Management

Name	Action
Configuration 1	Delete
Configuration 2	Delete
Configuration 3	Delete
Configuration 4	Delete
Configuration 5	Delete
Configuration 6	Delete
Configuration 7	Delete

[Create Configuration](#)

Please confirm.

[Yes](#)[No](#)

3

Modal for confirmation on delete.